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ROLL Long-Lasting, Drug-Free Pain Relief

PRODUCT DESCRIPTION

SportsWrap is a patented, hydrophilic, polyurethane membrane pad. The semi-permeable continuous thin film backing is optimized for moisture vapor permeability.

INDICATIONS

For the local management of swelling, inflammation, bruising, stiffness or pain due to acute or chronic injury when the skin is intact.

PRECAUTIONS

Consult your physician for all serious injuries, and if symptoms persist after wearing or increase during use. Not sterile. Not for broken skin or open wounds. For external use only.

PREPARATION

Clean the skin of dirt, topical medications and lotions. Residual moisture on skin is beneficial.

APPLICATION

Apply as soon as possible after injury. SportsWrap provides flexibility of application.

No matter which method of application you use or which part of the body you choose to wrap, the following steps will apply:

Step 1 Unroll the green pad with "THIS SIDE OUT" away from the skin.

Step 2 Wrap in a continuous motion overlapping slightly.

Step 3 Secure pad with athletic or adhesive tape.

Wrap firmly for close contact with the skin but apply gently to maintain comfort. The green pad must maintain intimate contact with the skin to work effectively.

TRIMMING

Depending on anatomy and wrapping technique, the SportsWrap may be too long. The green pad can be trimmed with a pair of scissors if necessary. Return the unused portion to the resealable pouch. Overlap to ensure complete coverage.

WEARING

1. Wear all day and night for 5 days.
2. Avoid wetting the wrap. Remove before bathing or showering.
3. If the green pad becomes saturated with sweat, discard and replace. The wrap cannot be "dried out" and reused.
4. Another SportsWrap may be worn for an additional 5 days. Continued use will help promote continued comfort.
5. Keep wearing even when you feel better!
6. Adjust periodically to maintain contact and comfort.

STORAGE

Keep dry and store at room temperature and humidity.

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Ferris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.

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© 2017, 2019 Ferris Mfg. Corp. ® or ™ indicate trademark